


Product Datasheet

DC-SIGN Monoclonal Antibody (orb75850)

Description	DC-SIGN Monoclonal Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Monoclonal
Uniprot ID	Q9NNX6
Expiration Date	12 months from date of receipt.



Western blot analysis of human uterus ti...



Immunohistochemical analysis of formalin...