


Product Datasheet

PUMA Monoclonal Antibody (orb75857)

Description	PUMA Monoclonal Antibody
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Monoclonal
Uniprot ID	Q96PG8
Expiration Date	12 months from date of receipt.



Western blot analysis of K562 cell lysat...