


## Product Datasheet

### Hemagglutinin Monoclonal Antibody (orb75865)

|                            |                                                                                           |
|----------------------------|-------------------------------------------------------------------------------------------|
| <b>Description</b>         | Hemagglutinin Monoclonal Antibody                                                         |
| <b>Species/Host</b>        | Mouse                                                                                     |
| <b>Reactivity</b>          | Virus                                                                                     |
| <b>Conjugation</b>         | Unconjugated                                                                              |
| <b>Tested Applications</b> | ELISA, WB                                                                                 |
| <b>Storage</b>             | Store at -20°C for long term preservation (recommended). Store at 4°C for short term use. |
| <b>Note</b>                | For research use only                                                                     |
| <b>Clonality</b>           | Monoclonal                                                                                |
| <b>Uniprot ID</b>          | <b>Q692M2</b>                                                                             |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                           |



Western blot analysis of 25 ng of recomb...