


Product Datasheet

Hax1a Monoclonal Antibody (orb75868)

Description	Hax1a Monoclonal Antibody
Species/Host	Mouse
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Monoclonal
Uniprot ID	000165
Expiration Date	12 months from date of receipt.



Western blot analysis of human brain tis...