



Product Datasheet

ORAI3 Monoclonal Antibody (orb75874)

Description	ORAI3 Monoclonal Antibody
Species/Host	Mouse
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Monoclonal
Uniprot ID	Q9BRQ5
Expiration Date	12 months from date of receipt.



Western blot analysis of rat spleen lysa...



Immunohistochemical analysis of formalin...