


Product Datasheet

XBP Monoclonal Antibody (orb75876)

Description	XBP Monoclonal Antibody
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Monoclonal
Uniprot ID	P17861
Expiration Date	12 months from date of receipt.



Western blot analysis of 100 ng of XBP-1...



Immunohistochemical analysis of formalin...