

Product Datasheet

IL-6R antibody (orb758841)

Description

Rabbit monoclonal antibody to IL-6R

Species/Host	Rat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	Blocking, FC
Immunogen	This antibody was raised by immunising Lewis rats with OKT4 cells.
Target	IL-6R
Preservatives	PBS with 0.02% Proclin 300.
Concentration	1 mg/ml
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Application notes	This antibody binds to IL-6R-positive cells, and blocks IL-6 binding (Vink et al, 1990; Coulie et al, 1990). This antibody is able to inhibit the growth of IL6-dependent plasmacytoma T1033C2 in vitro (Coulie et al, 1990). Administration of this antibody to mice prior to challenge with T1033C2 resulted in a protection rate of 60% and prolonged survival in tumour-bearing mice (Vink et al, 1990). This antibody has been used in flow cytometry to assess IL-6Ra expression (Hoge et al, 2013). When added to dendritic cells subject to CD40 activation in vitro, this antibody neutralises IL-6 and ablates the effect of CD40 activation (Grohmann et al, 2001).
Isotype	IgG
Clonality	Monoclonal
Purity	Purified
Uniprot ID	P22272
Expiration Date	12 months from date of receipt.