

Product Datasheet

IL-9 antibody (orb758847)

Description	Rabbit monoclonal antibody to IL-9
Species/Host	Rat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, FC
Immunogen	This antibody was raised by immunising rats with recombinant full-length murine IL-9.
Target	IL-9
Preservatives	PBS with 0.02% Proclin 300.
Concentration	1 mg/ml
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Application notes	This antibody, when fluorochrome-conjugated, has been used to achieve intracellular IL-9 staining of T cells, analysed using flow cytometry (Chang et al, 2010; Kim et al, 2010; Horka et al, 2012; Gomez-Rodriguez et al, 2016). This antibody has also been used in a sandwich ELISA to measure IL-9 secretion in supernatants collected from restimulated T cells (Jager et al, 2009).
Isotype	IgG
Clonality	Monoclonal
Purity	Purified
Uniprot ID	P15247
Expiration Date	12 months from date of receipt.