



Product Datasheet

Seasonal H1N1 Hemagglutinin Monoclonal Antibody (orb75888)

Description	Seasonal H1N1 Hemagglutinin Monoclonal Antibody
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Monoclonal
Uniprot ID	B0VX44
Expiration Date	12 months from date of receipt.



Western blot analysis of Recombinant seasonal influenza A H1N1 HA protein.



Seasonal Influenza A Hemagglutinin antibody specificity ELISA assay.