

Product Datasheet

Lewis a,b antibody (orb758887)

Description	Rabbit monoclonal antibody to Lewis a,b
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, IP, WB
Immunogen	This antibody was raised by immunising mice with human carcinoma cells.
Target	Lewis a,b
Preservatives	PBS with 0.02% Proclin 300.
Concentration	1 mg/ml
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Application notes	This antibody has applications in the detection of normal and neoplastic epithelial cells. This antibody has been shown to react with a glycolipid fraction isolated from colon carcinoma tissue, with no effect of pronase E addition on its reactivity (Moldenhauer et al, 1987). This antibody, when immobilised on the wall of a shear chamber, has been shown to support slow leukocyte rolling (Chen et al, 1997).
Isotype	IgG
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.