


Product Datasheet

Swine H1N1 Hemagglutinin Monoclonal Antibody (orb75889)

Description	Swine H1N1 Hemagglutinin Monoclonal Antibody
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Monoclonal
Uniprot ID	C4AL34
Expiration Date	12 months from date of receipt.



S-OIV A H1 Antibody specifically recogni...