

Product Datasheet

Human Amyloid Precursor protein (orb759025)

Description	Recombinant human Amyloid Precursor protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Lyophilized from a 0.2 um filtered solution of PBS, 5% trehalose, pH 7.4
Form/Appearance	Lyophilized from a 0.2 um filtered solution of PBS,pH7.4.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-Fc
Note	For research use only
Protein Sequence	Leu18-Lys612
Source	Human Cells
Uniprot ID	P05067
Expiration Date	6 months from date of receipt.