

Product Datasheet

Human GAS6 protein (orb759047)

Description	Recombinant human GAS6 protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 um filtered solution of PBS, 10% Glycerol , pH 7.4
Form/Appearance	Supplied as a 0.2 um filtered solution of PBS, 10% Glycerol , pH 7.4.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	C-6His
Note	For research use only
Protein Sequence	Ala31-Ala678
Source	Human Cells
Uniprot ID	Q14393
Expiration Date	6 months from date of receipt.