

Product Datasheet

Human IL36a protein (orb759099)

Description	Recombinant human IL36a protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Lyophilized from a 0.2 um filtered solution of 50mM Tris HCl, 150mM NaCl, pH 8.0,1 mM EDTA□0.02%Tween 20
Form/Appearance	Lyophilized from a 0.2 um filtered solution of 50mM Tris-HCl, 150mM NaCl, pH 8.0,1 mM EDTA,0.02%Tween 20.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	tag free
Note	For research use only
Protein Sequence	Lys6-Phe158
Source	E.coli
Uniprot ID	Q9UHA7
Expiration Date	6 months from date of receipt.