

## Product Datasheet

### Human CD73 protein (orb759147)

|                         |                                                                                                       |
|-------------------------|-------------------------------------------------------------------------------------------------------|
| <b>Description</b>      | Recombinant human CD73 protein                                                                        |
| <b>Endotoxins</b>       | Less than 0.1 ng/ug (1 EU/ug)                                                                         |
| <b>Preservatives</b>    | Supplied as a 0.2 um filtered solution of 20mM TrisHCl, 120mM NaCl, 4mM CaCl2, 20% Glycerol, pH 7.5   |
| <b>Form/Appearance</b>  | Supplied as a 0.2 um filtered solution of 20mM Tris-HCl, 120mM NaCl, 4mM CaCl2, 20% Glycerol, pH 7.5. |
| <b>Storage</b>          | Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.                |
| <b>Tag</b>              | C-Fc                                                                                                  |
| <b>Note</b>             | For research use only                                                                                 |
| <b>Protein Sequence</b> | Trp27-Lys547                                                                                          |
| <b>Source</b>           | Human Cells                                                                                           |
| <b>Uniprot ID</b>       | AAH65937.1                                                                                            |
| <b>Expiration Date</b>  | 6 months from date of receipt.                                                                        |