

Product Datasheet

Human CD73 protein (orb759165)

Description	Recombinant human CD73 protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 um filtered solution of 20mMTrisHCl,120mMNaCl,4mM CaCl_2 ,20%Glycerol,pH7.5
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mMTris-HCl,120mMNaCl,4mM CaCl_2 ,20%Glycerol,pH7.5.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	C-6His
Note	For research use only
Protein Sequence	Trp27-Lys547
Source	Human Cells
Uniprot ID	P21589
Expiration Date	6 months from date of receipt.