

Product Datasheet

Primate CD16a protein (orb759169)

Description	Recombinant primate CD16a protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Lyophilized from a 0.2 um filtered solution of PBS, pH 7.4
Form/Appearance	Lyophilized from a 0.2 um filtered solution of PBS, pH 7.4.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-6His
Note	For research use only
Protein Sequence	Glu21-Gly206
Source	Human Cells
Uniprot ID	A0A140HDP8
Expiration Date	6 months from date of receipt.