

Product Datasheet

Human IL23 protein (orb759171)

Description	Recombinant human IL23 protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Lyophilized from a 0.2 um filtered solution of 20mMPB,150mMNaCl,pH7.4
Form/Appearance	Lyophilized from a 0.2 um filtered solution of 20mMPB,150mMNaCl,pH7.4.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-Fc
Note	For research use only
Protein Sequence	Ile23-Ser328&Ala21-Pro189
Source	Mammalian
Uniprot ID	Q9NPF7, P29460
Expiration Date	6 months from date of receipt.