

Product Datasheet

Human CD30 protein (Biotin) (orb759180)

Description	Recombinant human CD30 protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Conjugation	Biotin
Preservatives	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mMNaCl, pH 7.4.
Form/Appearance	Lyophilized from a 0.2 um filtered solution of 20mM PB, 150mMNaCl, pH 7.4.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-6His-Avi
Note	For research use only
Protein Sequence	Gln24-Lys265(Met115Thr)
Source	Human Cells
Uniprot ID	P28908
Expiration Date	6 months from date of receipt.