

Product Datasheet

Human IL-2 R gamma / CD132 Protein (orb759218)

Description

Human IL-2 R gamma, His Tag is

(HEK293). It contains AA Leu 23 - Asn ...

Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

IL-2 R gamma

Preservatives

PBS, pH7.4

Form/Appearance

Powder

Storage

-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Tag

C-10×His

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 29.3 kDa. The protein migrates as 55-65 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Purity

95%

Source

Human IL-2 R gamma, His Tag (orb759218) is expressed from human 293 cells (HEK293). It contains AA Leu 23 - Asn 254 (Accession # P31785-1).

MW

29.3 kDa

Uniprot ID

P31785

NCBI

AAH14972.1

Expiration Date

6 months from date of receipt.