


Product Datasheet

Free Fatty Acid Microplate Assay Kit (orb759221)

Description	Detection and Quantification of Free Fatty Acid Content.
Reactivity	All
Range	0.01 mmol/L - 1 mmol/L
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Serum, Plasma, Urine, Saliva, Milk, Cell culture medias, Other biological fluids, Platform: Microplate reader 550 nm
Sample Types	Serum, Plasma, Urine, Saliva, Milk, Cell culture medias, Other biological fluids
Expiration Date	Please enquire.



Representative
Standard
Curve