


Product Datasheet

Total Bile Acid Microplate Assay Kit (orb759222)

Description	Detection and Quantification of Bile Acid Content.
Reactivity	All
Range	2 µmol/L - 200 µmol/L
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Urine, Serum, Tissue extracts, Other biological fluids, Platform: Microplate reader 540 nm
Sample Types	Urine, Serum, Tissue extracts, Other biological fluids
Expiration Date	Please enquire.



Representative
Standard
Curve