


Product Datasheet

Histamine Microplate Assay Kit (orb759223)

Description	Detection and Quantification of Histamine Content.
Reactivity	All
Range	0.02 mmol/L - 2 mmol/L
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Tissue extracts, Other biological fluids, Platform: Microplate reader 480 nm
Sample Types	Tissue extracts, Other biological fluids
CAS Number	51-45-6
Expiration Date	Please enquire.



Representative
Standard
Curve