

Product Datasheet

Tannase Microplate Assay Kit (orb759225)

Description Detection and Quantification of Tannase Activity.

Reactivity All

Range 0.02 mmol/L - 2 mmol/L

Tested Applications FA

Storage Shipped and store at 4°C for 6 months.


Note For research use only

Application notes Detection Method: Colorimetric, Sample Types: Tissue extracts, Cell lysate, Cell culture media, Other biological fluids, Platform: Microplate reader 520 nm


Sample Types Tissue extracts, Cell lysate, Cell culture media, Other biological fluids

CAS Number 9025-71-2

Expiration Date Please enquire.



Representative
Standard
Curve



Positive
Control
reaction in 96-
well pla...