


Product Datasheet

Indoleacetic Acid Oxidase Microplate Assay Kit (orb759226)

Description	Detection and Quantification of Indoleacetic Acid Oxidase Activity.
Reactivity	All
Range	0.1 mmol/L - 2 mmol/L
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Tissue extracts, Cell lysate, Cell culture media, Other biological fluids, Platform: Microplate reader 530 nm
Sample Types	Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
Expiration Date	Please enquire.



Representative
Standard
Curve