


Product Datasheet

Galactose Microplate Assay Kit (orb759227)

Description	Detection and Quantification of Galactose Content.
Reactivity	All
Range	0.01 mmol/L - 1 mmol/L
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Serum, Plasma, Urine, Saliva, Milk, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids, Platform: Microplate reader 450 nm
Sample Types	Serum, Plasma, Urine, Saliva, Milk, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
CAS Number	59-23-4
Expiration Date	Please enquire.



Representative
Standard
Curve