

Product Datasheet

Lactose Microplate Assay Kit (orb759228)

Description Detection and Quantification of Lactose Content.

Reactivity All

Range 0.2 mmol/L - 20 mmol/L

Tested Applications FA

Storage Shipped and store at 4°C for 6 months.


Note For research use only

Application notes Detection Method: Colorimetric, Sample Types: Serum, Plasma, Urine, Saliva, Milk, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids, Platform: Microplate reader 505 nm

Sample Types Serum, Plasma, Urine, Saliva, Milk, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids

CAS Number 63-42-3

Expiration Date Please enquire.



Representative
Standard
Curve