


Product Datasheet

Pyruvate Carboxylase Microplate Assay Kit (orb759229)

Description	Detection and Quantification of Pyruvate Carboxylase Activity.
Reactivity	All
Range	4 µmol/L - 400 µmol/L
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Urine, Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids, Platform: Microplate reader 340 nm
Sample Types	Urine, Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
CAS Number	9014-19-1
Expiration Date	Please enquire.



Representative
Standard
Curve