



Product Datasheet

Uricase Microplate Assay Kit (orb759230)

Description	Detection and Quantification of Uricase Activity.
Reactivity	All
Range	0.03 mmol/L - 3 mmol/L
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Tissue extracts, Cell lysate, Cell culture media, Other biological fluids, Platform: Microplate reader 550 nm
Sample Types	Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
CAS Number	9002-12-4
Expiration Date	Please enquire.



Representative Standard Curve



Positive Control reaction in 96-well pla...