

## Product Datasheet

### Broad Spectrum Phosphatase Inhibitor Cocktail [Out of Stock] [Out of stock] (orb76033)

<b>Description</b>	Broad Spectrum Phosphatase Inhibitor Cocktail is a Western blot related concentrated stock solution ...
<b>Tested Applications</b>	WB
<b>Storage</b>	Upon receipt store at -20°C. It is stable at -20°C for one year. It is stable at 4°C for two months. Product is shipped on ice.
<b>Note</b>	For research use only
<b>Application notes</b>	WB sample preparation, protein purification, Co-IP, protein kinase activity assay experiments
<b>Expiration Date</b>	12 months from date of receipt.