


Product Datasheet

Human CCL24/Eotaxin-2 ELISA Kit (orb76036)

Description	Human CCL24/Eotaxin-2 ELISA Kit
Reactivity	Human
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: V27-C119
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	13134 MW
Uniprot ID	000175
Sensitivity	2 pg/ml
Dilution Range	Dilution ratio of 1:5, concentration in serum and plasma is 188-3727 pg/ml.
Expiration Date	Please enquire.



Human CCL24/Eotaxin-2 PicoKine ELISA Kit...