


Product Datasheet

Human FASL / CD95LG / CD95 Ligand ELISA Kit (orb76037)

Description	Human FASL / CD95LG / CD95 Ligand ELISA Kit
Reactivity	Human
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: CHO; Immunogen sequence: P134-L281
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA, citrate).
MW	31485 MW
Uniprot ID	P48023
Sensitivity	2 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is 40-218 pg/ml.
Expiration Date	Please enquire.



Human
FASL
PicoKine
ELISA Kit
standard
C...