


Product Datasheet

Human CD25/IL-2sR Alpha ELISA Kit (orb76042)

Description	Human CD25/IL-2sR Alpha ELISA Kit
Reactivity	Human
Range	78 pg/ml - 5,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: E22-C213
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	30819 MW
Uniprot ID	P01589
Sensitivity	5 pg/ml
Dilution Range	Dilution ratio of 1:50, concentration in serum and plasma is 670-1178 pg/ml.
Expiration Date	Please enquire.



Human
CD25/IL-2sR
alpha
PicoKine
ELISA K...