


Product Datasheet

Mouse Leptin Receptor ELISA Kit (orb76045)

Description	Mouse Leptin Receptor ELISA Kit
Reactivity	Mouse
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: L22-G839
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	130789 MW
Uniprot ID	P48356
Sensitivity	5 pg/ml
Expiration Date	Please enquire.



Mouse
Leptin
receptor
PicoKine
ELISA Kit...