


Product Datasheet

Human CCL7/MCP-3 ELISA Kit (orb76047)

Description	Human CCL7/MCP-3 ELISA Kit
Reactivity	Human
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: Q24-L99
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	11200 MW
Uniprot ID	P80098
Sensitivity	1 pg/ml
Expiration Date	Please enquire.



Human
CCL7/MCP-3
PicoKine
ELISA Kit
stan...