


Product Datasheet

Human VEGFR2 / KDR / VEGF Receptor 2 ELISA Kit (orb76050)

Description	Human VEGFR2 / KDR / VEGF Receptor 2 ELISA Kit
Reactivity	Human
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: A20-E764
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	151527 MW
Uniprot ID	P35968
Sensitivity	4 pg/ml
Dilution Range	Dilution ratio of 1:10, concentration in serum and plasma is 5136-17401 pg/ml.
Expiration Date	Please enquire.



Human
VEGFR2/KDR
PicoKine
ELISA Kit
stan...