


Product Datasheet

Human CCL4/MIP-1 Beta ELISA Kit (orb76052)

Description	Human CCL4/MIP-1 Beta ELISA Kit
Reactivity	Human
Range	7.8 pg/ml - 500 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: A24-N92
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA, citrate).
MW	10212 MW
Uniprot ID	P13236
Sensitivity	1 pg/ml
Dilution Range	Dilution ratio of 1:4, concentration in serum and plasma is 200-400 pg/ml.
Expiration Date	Please enquire.



Human
CCL4/MIP-1
beta
PicoKine
ELISA Kit...