


Product Datasheet

Human CCL21 / 6CKine / Exodus 2 ELISA Kit (orb76053)

Description	Human CCL21 / 6CKine / Exodus 2 ELISA Kit
Reactivity	Human
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: S24-P134
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	14646 MW
Uniprot ID	000585
Sensitivity	12 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is around 500 pg/ml.
Expiration Date	Please enquire.



Human
CCL21/6CKine
PicoKine
ELISA Kit st...