


Product Datasheet

Human C-MET/HGFR ELISA Kit (orb76058)

Description	Human C-MET/HGFR ELISA Kit
Reactivity	Human
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: E25-R307 (alpha) & S308-T932 (beta)
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	155541 MW
Uniprot ID	P08581
Sensitivity	5 pg/ml
Expiration Date	Please enquire.



Human C-MET/HGFR PicoKine ELISA Kit stan...