


Product Datasheet

Human KLK5/Kallikrein-5 ELISA Kit (orb76061)

Description	Human KLK5/Kallikrein-5 ELISA Kit
Reactivity	Human
Range	78 pg/ml - 5,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: V23-S293
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA) , human milk and saliva.
MW	32019 MW
Uniprot ID	Q9Y337
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is around 140 pg/ml.
Expiration Date	Please enquire.



Human
KLK5
PicoKine
ELISA Kit
standard c...