


Product Datasheet

Human TLR3/Toll-like receptor 3 ELISA Kit (orb76064)

Description	Human TLR3/Toll-like receptor 3 ELISA Kit
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: K27-S711
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	103829 MW
Uniprot ID	O15455
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is less than the lowest standard, 156 pg/ml.
Expiration Date	Please enquire.



Human
TLR3
PicoKine
ELISA Kit
standard
C...