

Product Datasheet

Human COMP/Thrombospondin 5 ELISA Kit (orb76065)

Description Human COMP/Thrombospondin 5 ELISA Kit

Reactivity Human

Range 156 pg/ml - 10,000 pg/ml

Tested Applications ELISA

Immunogen Expression system for standard: NS0; Immunogen sequence: Q21-A757

Storage Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)

Note For research use only

Sample Types cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).


MW 82860 MW

Uniprot ID **P49747**

Sensitivity 10 pg/ml

Dilution Range Dilution ratio of 1:50, concentration in serum and plasma is 53.44-490.8 ng/ml.

Expiration Date Please enquire.



Human
COMP
PicoKine
ELISA Kit
standard c...