

Product Datasheet

Lamin A antibody (FITC) (orb761194)

Description	Mouse monoclonal antibody against Lamin A conjugated to FITC
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	FITC
Tested Applications	IF, IHC-Fr
Immunogen	Recombinant human lamin A
Target	Lamin A
Preservatives	Purified IgG conjugated to Fluorescein using Innova Biosciences LightningLink, supplied in Phosphate buffered saline (PBS) containing 0.09% Sodium azide
Concentration	1mg/ml
Note	For research use only
Application notes	Spleen cells from immunised BALB/c mice were fused with cells of the mouse SP2-0/Ag4 myeloma cell line
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P02545
Entrez	4000
Dilution Range	Optimal antibody dilution should be determined by titration
Expiration Date	12 months from date of receipt.