

Product Datasheet

Somatostatin Receptor 1/SSTR1 Antibody (orb76207)

Description

Somatostatin Receptor 1/SSTR1 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

IHC, WB

Applications

Immunogen

A synthetic peptide corresponding to a sequence at the C-terminus of human Somatostatin Receptor 1 (361-381aa RAYSEDFQPENLESGGVFRN), identical to the related rat and mouse sequences.

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Predicted Species: Species predicted to be fit for the product based on sequence similarities. By Heat: Boiling the paraffin sections in 10mM citrate buffer, pH6.0, for 20mins is required for the staining of formalin/paraffin sections. Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

42686 MW

Uniprot ID

P30872

Dilution Range

Immunohistochemistry (Paraffin-embedded Section), 0.5-1µg/ml, Rat, Human, Mouse, By Heat Western blot, 0.1-0.5µg/ml, Rat, Human, Mouse

Expiration Date

12 months from date of receipt.

100KD —


55KD —

35KD —

25KD —

15KD —

Anti-SSTR1 antibody, orb76207, Western b...



Anti-SSTR1 antibody, orb76207, IHC (P) I...