

Product Datasheet

Cynomolgus MCAM / CD146 Protein (orb762441)

Description	Cynomolgus MCAM, His Tag is expressed from human 293 cells (HEK293). It contains AA Val 24 - Gly 559...
Reactivity	Monkey
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	MCAM
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 61.5 kDa. The protein migrates as 80-110 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	XP_015291315.1
Purity	95%
Source	Cynomolgus MCAM, His Tag (orb762441) is expressed from human 293 cells (HEK293). It contains AA Val 24 - Gly 559 (Accession # A0A5F8A4A1-1).
MW	61.5 kDa
Uniprot ID	A0A5F8A4A1
NCBI	XP_015291315.1
Expiration Date	6 months from date of receipt.