


Product Datasheet

Human IFNA2 ELISA Kit (orb763006)

Description	Human IFNA2 ELISA Kit
Reactivity	Human
Range	3.12 pg/ml - 200 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: C24-E189
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants and serum
Uniprot ID	P01563
Sensitivity	1 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum is less than the lowest standard, 3.1 pg/ml.
Expiration Date	Please enquire.



Human
IFNA2
PicoKine
ELISA Kit
standard ...