


Product Datasheet

Human Pro-Cathepsin B ELISA Kit (orb763008)

Description	Human Pro-Cathepsin B ELISA Kit
Reactivity	Human
Range	0.156 ng/ml - 10 ng/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: R18-I339
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin), saliva and urine
Uniprot ID	P07858
Sensitivity	20 pg/ml
Expiration Date	Please enquire.



Human Pro-Cathepsin B PicoKine ELISA Kit...