

Product Datasheet

Grp75/HSPA9 Antibody (orb76302)

Description

Grp75/HSPA9 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested Applications

ICC, IHC, IHC-Fr, WB

Immunogen

A synthetic peptide corresponding to a sequence at the C-terminus of human Grp75 (591-609aa DTETKMEEFKDQLPADECN), identical to the related rat and mouse sequences.

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

73680 MW

Uniprot ID

P38646


Dilution Range

Western blot, 0.1-0.5µg/ml, Human, Mouse, Rat
Immunohistochemistry (Paraffin-embedded Section), 0.5-1µg/ml, Human, Mouse, Rat, By Heat
Immunohistochemistry (Frozen Section), 0.5-1µg/ml, Rat Immunocytochemistry , 0.5-1µg/ml, Human


Expiration Date

12 months from date of receipt.

1 2 3 4 5



Western blot analysis of Lane 1: Rat Liv...



Immunohistochemical analysis of formalin...

