

Product Datasheet

CAMKIV/Camk4 Antibody (orb763026)

Description

CAMKIV/Camk4 Antibody

Species/Host

Rabbit

Reactivity

Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, FC, ICC, IF, IHC, WB

Applications
Immunogen

E.coli-derived mouse Camk4 recombinant protein (Position: A19-Y469).

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW


58 kDa

Uniprot ID
P08414
Dilution Range


Western blot, 0.1-0.25µg/ml, Mouse, Rat
Immunohistochemistry (Paraffin-embedded Section), 2-5µg/ml, Mouse, Rat
Immunocytochemistry/Immunofluorescence, 5µg/ml, Mouse
Immunofluorescence, 5µg/ml, Rat Flow Cytometry, 1-3µg/1x10⁶ cells, Mouse Direct ELISA, 0.1-0.5µg/ml, Mouse

Expiration Date


12 months from date of receipt.



Western blot analysis of CAMKIV/Camk4 using anti-CAMKIV/Camk4 antibody.



IHC analysis of CAMKIV/Camk4 using anti-CAMKIV/Camk4 antibody.



IHC analysis of CAMKIV/Camk4 using anti-CAMKIV/Camk4 antibody.