



Product Datasheet

Ki67/Mki67 Antibody (orb763060)


Description	Ki67/Mki67 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	E.coli-derived mouse Mki67 recombinant protein (Position: F29-S3177).
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Add 0.2ml of distilled water will yield a concentration of 500ug/ml.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	358 kDa
Uniprot ID	E9PVX6
Dilution Range	Immunohistochemistry (Paraffin-embedded Section), 2-5µg/ml, Human, Mouse, Rat Direct ELISA, 0.1-0.5µg/ml, Mouse
Expiration Date	12 months from date of receipt.



IHC analysis of Ki67/Mki67 using anti-Ki...



IHC analysis of Ki67/Mki67 using anti-Ki...



IHC analysis of Ki67/Mki67 using anti-Ki...